
NL405 - Width/Delay



Introduction

The **NL405 WIDTH/DELAY** is a digitally controlled pulse duration generator and is used to produce precise pulses of adjustable duration with high accuracy and repeatability.

The **NL405** produces TTL logic output pulses (of both polarities) that have a duration set using three decades of thumbwheels and a seven (7) position range switch. It also features a press to start button, LED indication of an output and a separate logic output signal at the end of the defined WIDTH time. It was designed as a replacement for the NL401 and in most places where the NL403 would have been used.

The **NL405** can be triggered off either positive or negative-going edges of the input pulses and the **NL405** is capable of 100% duty cycles. An external push to make switch can be used when -ve edge triggering is selected.

In the "FIXED WIDTH" mode, the **NL405** produces an Output pulse as specified by the front panel controls; In the "RE-TRIGGER" mode, any valid Input trigger arriving whilst the Output is active will restart the timing and therefore extending the Output pulse duration.

As long time periods can be generated (999 seconds is over 16 minutes) any change in position of the rotary range switch will stop the current Output and inhibit further input triggers for 450 ms.

The "DELAYED" Output socket provides a short duration pulse (set internally) that is produced at the end of the timed WIDTH pulse. This provides a trigger-out after the set Delay time.


The **NL405** can be used for digital control of pulse duration or for introduction of a delay to a trigger pulse. The logarithmic or exponential nature of the many biological responses make the very wide continuous adjustment (from μ s to minutes) particularly beneficial.

Specification Summary

Input:	:	TTL pulse. Input is 1 TTL load.
Trigger edge	:	Positive or Negative, selected from front panel
Manual	:	Push-button, operation on push.
Output (Width)	:	TTL. Output will drive 10 loads.
Polarity	:	Both active-high and active-low sockets
Range	:	1 μ s to 999 s
Control	:	001 to 999
Multiplier	:	1 μ s, 10 μ s, 100 μ s, 1 ms, 10 ms, 100 ms, 1 s.
Accuracy	:	$\pm 0.5\%$
Delay from Trigger In	:	<200 ns
Indicator	:	Amber LED on for approximately 100 ms for each pulse
Output (Delay)	:	TTL. Output will drive 10 loads.
Polarity	:	Active-high
Duration	:	0.5 μ s, 1.5 μ s, 5 μ s $\pm 20\%$ 50 μ s, 150 μ s, 500 μ s $\pm 10\%$, by on-board jumpers.
Output		
Mode	:	Fixed Width or Re-Triggered, selected from front panel
Reset/Inhibit	:	450 ms $\pm 20\%$, after change of rotary range control.

Rear connections to the motherboard allow Input and Output interconnections between this and other logic modules without the need of front panel cables.


Mode selection





In **FIXED WIDTH** mode, the NL405 will produce an Output pulse of the duration as specified on the Front Panel regardless of further Trigger applied to the input (IN). Trigger pulses received during the time the Output (OUT) is active are ignored.

In **RE-TRIGGER** mode, the Output pulse duration of NL405 is extended when a further trigger pulse is received at the input (IN). The number of times an output can be extended by re-triggering, is unlimited.

Output "reset" and Inhibit



The Output pulse duration can be terminated early, in either mode, by rotating the Range control. It would be normal, in this instance, to rotate it to the next 'click' position and then back to where it was set. The Outputs are inhibited for 450 ms after the switch becomes still. Note that a Delayed Output pulse is NOT generated at the end of a WIDTH that is terminated early.


Trigger selection and relationship to Outputs

The above two figures show the effect of selecting “+ve edge” or “-ve edge” triggering.

POSITIVE EDGE (+ve edge) triggering means that the unit will be internally triggered when a rising edge is seen on the input (IN) socket.

NEGATIVE EDGE (-ve edge) triggering means that the unit will be internally triggered when a falling edge is seen on the input (IN) socket.

This can often be used with an external switch closure, such as a foot switch, as long as the DELAYED OUT is used, the unit is set to RE-TRIGGER and the WIDTH is set longer than the contact bounce time (which will normally be under 10 ms). The unit may also be triggered when the switch is released (because of poor contacts) and in this instance an external capacitor across the contacts is called for. This may need to be found from trial and error but a value of 1 μ F would be a good starting point.

Rear Connection IN/OUT jumpers (on-board)

The rear edge connector in the NL900 rack allows adjacent modules to be connected together without the need of external front panel leads.


This module has a choice of three output signals that can be connected to the rear connector for automatic routing to the module on the immediate right. The Input signal can be disconnected from the output of the module on the immediate left if it is inappropriate.

The position of these jumpers are changed by gently pulling them away from the PCB, then rotating them 90°, 180° or 270° before replacing them in two of the small gold board sockets. In the "Parked" or "N/C" position the jumper is placed in only one gold socket and a parking hole.

INPUT - the upper jumper allows the Output signal from the module to the immediate left to be connected to the input socket (IN), or it may be "parked" for "no-connection (N/C)".

OUTPUT - the lower jumper selects which of the three Output signals is passed to the Input of the module to the immediate right, or it may be "parked" for "no-connection (N/C)".

Delayed Output Pulse Duration Jumpers (on-board)



Delayed Output Pulse Duration Selection

Jumpers	Setting	Jumpers	Setting
	0.5 μs		50 μs
	1.5 μs		150 μs
	5 μs		500 μs

The Duration of the pulse at “DELAYED OUT” can be set to one of six (6) pulse durations and is set by two jumpers on the PCB.

These are altered by pulling them off their pins, moving them to where is required and pushing them back on. Both jumpers **MUST** be fitted in one of the six configurations shown below.

We reserve the right to alter specifications and price without prior notification.

First Issued: December 19, 2002

Last Revision: March 26, 2003

Printed: March 26, 2003

File Reference: N:\Docs\Company\Manuals\NeuroLog\NL405-02.lwp