



Matlab Control of the Digitimer DS8R Stimulator

Introduction

The Digitimer DS8R Stimulator incorporates a Windows compatible software API which enables third party control of the stimulator via Matlab, Python and other software packages. Here we provide basic instructions for implementing Matlab control and include a sample program that provides a basic demonstration of this capability.

Hardware & Software Requirements

- DS8R Stimulator.
- USB cable (supplied with DS8R).
- Personal Computer running 32bit or 64bit version of Windows 10.
- Digitimer DS8R Virtual Front Panel Software Installer (supplied with DS8R and available [here](#)).
- Digitimer DS8R Matlab Solution ZIP file, including the DS8R Proxy DLL installer, available [here](#).
- Base Matlab package, including a third-party compiler.

Installation Process

1. Download, install and run the DS8R Virtual Front Panel Software.
2. Connect the DS8R to the PC via the USB cable. Power ON the DS8R and verify that the DS8R Virtual Front Panel Software is able to control the DS8R settings.
3. Download and extract the Digitimer Matlab Solution ZIP file.
4. Locate and run D128RProxyv2.exe from within the DS8R Proxy DLL Installer subfolder of the ZIP file contents. This installs the DS8R Proxy DLL (D128RProxy.dll) into the correct system folder (C:\Windows\System32 or C:\Windows\SysWOW64) for Windows 32bit and Windows 64bit installations.
5. Copy the remaining three files (D128ctrl.m, D128RProxy.h, Ex_D128.m) to your chosen working directory for Matlab projects.

Demonstrating Matlab Control

Run the Matlab application and open the working folder containing the three downloaded files:-

D128ctrl.m & D128RProxy.h

Together these two files create a Matlab interface to Digitimer DS8R devices which provides an easy-

to-use method for loading and calling functions in the previously installed DS8R Proxy DLL.

Ex_D128.m

This simple example program demonstrates control of the DS8R.

```
% Simple example program to demonstrate control of D128

% open device and return handle for further calss
[success, d128] = D128ctrl('open');

% Download status from device
[success, d128] = D128ctrl('status', d128);

[success, d128] = D128ctrl('enable', d128, 0);

% Set value of pulsewidth, but does not upload to device
[success, d128] = D128ctrl('source', d128, 'Internal');
[success, d128] = D128ctrl('pulsewidth', d128, 1000);
[success, d128] = D128ctrl('demand', d128, 60);
[success, d128] = D128ctrl('dwell', d128, 400);

% Uploads all parameters to device
success = D128ctrl('upload', d128);

[success, d128] = D128ctrl('enable', d128, 1);

% trigger the device
success = D128ctrl('Trigger', d128);

% Download status from device
[success, d128] = D128ctrl('status', d128);
d128

% Close device
success = D128ctrl('close', d128);
```

By running the “ex_d128.m” sample program, the settings below are applied to the DS8R.

- Pulse Duration 1000 μ s
- Current Amplitude 60mA
- Interphase Interval 500 μ s

The program then enables the DS8R output and triggers delivery of a single stimulus pulse.

Acknowledgements

We gratefully acknowledge the assistance provided by John Kakareka, Troy Dildine & Qingbao Yu (National Institutes of Health, USA), who collaborated with Digitimer over the development of this Matlab control interface.

Digitimer Limited

37 Hydeway, Welwyn Garden City, AL7 3BE, UK
Tel: +44(0)1707 328347 Email: sales@digitimer.com

Last Updated 5/10/21