

D177 Bio-Feedback Unit

More Reproducible VEMP Testing • Portable & Lightweight • Easy to Setup & Use



Designed to Improve the Reproducibility of Neurophysiological Testing

The Digitimer D177 Bio-Feedback Unit has been developed to improve reproducibility of the Vestibular Evoked Myogenic Potential (VEMP) test, however, it may be employed for any neurophysiological evaluation that requires repetitive voluntary input from the subject. The D177 provides a visual indication of the effort being exerted and this helps the subject maintain the same level of tonic activation during each successive trial.

Easily Incorporated into Standard EP/EMG Systems.

The D177 has a 2.8m long cable, terminating with a 3.5mm mono audio jack, to facilitate

connection to the audio output of an EMG amplifier. The D177 converts this signal into visual feedback for the subject, in the form of needle movement on an analogue level meter. The subject is instructed to make an initial movement under the direction of the operator and the 'volume' is then adjusted until the needle is in the centre of the display. No further adjustments are necessary until the test is complete for that person. All the operator need do, is ask the subject to exert enough voluntary effort to direct the needle to the centre of the display and maintain it there, during each trial.

Inexpensive and Suitable for Subjects with Large or Small EMG Responses.

The D177 Biofeedback Unit is a simple, inexpensive device which improves the reproducibility of neurophysiological tests including VEMPs. It is suitable for subjects with small or large EMG output, as the device is easily calibrated to an effort that each person can withstand and repeat comfortably.

For more information about the D177 Bio-Feedback Unit, please contact:-

Digitimer Ltd

37 Hydeway, Welwyn Garden City, Hertfordshire, AL7 3BE, UK

Tel: +44 [0]1707 328347 Fax: +44 [0]1707 373153

E-mail: sales@digitimer.com

Website: www.digitimer.com