

## D360 Audio Interface (Switched)

---



The D360-AIS gives the Operator a tool to select which channel to listen to or feedback to the patient.

The D360-AIS plugs into the Signal Channel Output socket on the rear of the D360 and the LINE-IN socket of a PC.



The PC (assuming a good audio card) is used to control the frequency response (tone - Bass & Treble) as well as volume.

The D360-AIS also has a volume control.

The D360-AIS contains audio-isolation transformers to remove the possibility of mains-frequency "hum" on the monitored signal.

The D360-AIS is supplied (as shown above) with an attached lead to the D360 Signal Output socket as well as a 3.5mm stereo lead to the PC "LINE-IN" socket (often coloured - Blue).

---

**Digitimer Ltd.** 37 Hydeway, Welwyn Garden City, Hertfordshire, AL7 3BE. UK

Tel: +44 (0)1707 328347 • Fax: +44 (0)1707 373153 • sales@digitimer.com • www.digitimer.com